



## Automatic Manosonic Nebulizing System

### Treatment of tubo-tympanic pathologies:

- secretory otitis media or seromucous otitis,
- tubotympanic sequelae (like retraction pocket...).

3 combined functions, automatically synchronized with the patient's swallowing, allow the aerosol penetration into the Eustachian Tube:

- aerosol production,
- emission of 100 Hz sonic vibrations,
- temporary application of a modulable adjustable overpressure.

### Contra-indications:

- child under 4 years old,
- patient with a patulous Eustachian tube,
- patient with cleft lip and palate.



Latex free accessories: **PURENEB® EXTRANEB MANOSONIC nebulizing kits** with overpressure tight nose pieces.



**EXTRANEB MANOSONIC** adult, ref. MS1A  
medium and large nose pieces



**EXTRANEB MANOSONIC** paediatric, ref. MS1E  
small and medium nose pieces



**Carrying bag**, ref. F64



**ATOMISOR CONTROL<sup>1</sup>** to check the good  
functioning of your device, ref. MEA1

*Class IIa Medical Devices. Carefully read the instructions for use provided with the medical device. Documentation for Health Professionals. Non contractual pictures*



A. COMPS.EN.044.0120