

Ordering & additional information

Device specs

			
	ECHELON LINEAR™ Cutter	ETHICON™ Linear Cutter	PROXIMATE™ Linear Cutter
Product Codes	GLC60, GLC80, GLC100	NTLC55, NTLC75	TCD75, TCT10, TCT55, TCT75, TLC10, TLC55, TLC75, TVC55
Associated Reloads	GLCR60B (60mm, Blue) GLCR60G (60mm, Green) GLCR80B (80mm, Blue) GLCR80G (80mm, Green) GLCR100B (100mm, Blue) GLCR100G (100mm, Green)	SR55 SR75	TVR55, TCR55, TRT55, TCR75 TRD75, TRT75, TCR10, TRT10
Description	Linear Cutter with Reloads Sold Separately	Selectable Linear Cutter with 3D Staples	Linear Cutter with Reload
Device Technology	3D Stapling and Gripping Surface Technology	Selectable Staple Height	Cam Closure Mechanism
Size	60mm, 80mm, 100mm	55mm, 75mm	55mm, 75mm, 100mm
Staple Rows	4	6	4
Closed Staple Height	1.5mm, 2.0mm	1.5mm, 1.8mm, 2.0mm	1.5mm, 1.8mm, 2.0mm
Max Firings	8	12	8
Staple Technology	3D	3D	B-Form
Qty./Box	3 staplers per box 12 reloads per box	3 staplers per box 12 reloads per box	3 staplers per box 12 reloads per box